


Product datasheet for **TA384397M**

HARS1 Mouse Monoclonal Antibody [Clone ID: 4H5-D9-D4]

Product data:

Product Type:	Primary Antibodies
Clone Name:	4H5-D9-D4
Applications:	WB
Recommended Dilution:	WB: 1/2000
Reactivity:	Human, Mouse, Monkey
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Purified recombinant human HARS protein fragments expressed in E.coli.
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.03% Proclin 300, pH 7.3.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	56kDa
Gene Name:	histidyl-tRNA synthetase
Database Link:	Entrez Gene 3035 Human P12081
Background:	Swiss-Prot Acc.P12081.
Synonyms:	FLJ20491; HisRS; HRS

[View online »](#)

Product images:


Western blot detection of HARS in T47D,3T3,COS7,C6 and HeLa cell lysates using HARS mouse mAb(dilution 1:2000).Predicted band size:56kDa.Observed band size:56kDa.