

Product datasheet for **TA384399M**


TIM 1 (HAVCR1) Rabbit Monoclonal Antibody [Clone ID: R03-6E8]

Product data:


Product Type:	Primary Antibodies
Clone Name:	R03-6E8
Applications:	WB
Recommended Dilution:	WB: 1/2000-1/10000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A synthetic peptide of human TIM 1
Formulation:	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Calculated MW:39 kDa,Observed MW:50 kDa
Gene Name:	hepatitis A virus cellular receptor 1
Database Link:	Entrez Gene 26762 Human Q96D42
Background:	The protein encoded by this gene is a membrane receptor for both human hepatitis A virus (HHAV) and TIMD4. The encoded protein may be involved in the moderation of asthma and allergic diseases. The reference genome represents an allele that retains a MTTVP amino acid segment that confers protection against atopy in HHAV seropositive individuals. Alternative splicing of this gene results in multiple transcript variants. Related pseudogenes have been identified on chromosomes 4, 12 and 19. [provided by RefSeq, Apr 2015]
Synonyms:	HAVCR; HAVCR-1; KIM-1; KIM1; TIM; TIM-1; TIM1; TIMD-1; TIMD1


[View online »](#)

Product images:



Western blot analysis of TIM 1 in Raw264.7, PC-12 lysates using Hepatitis A Virus ular Receptor 1 antibody.



Western blot analysis of TIM 1 in A549, MCF-7, U251, HL-60 lysates using TIM 1 antibody.