

## Product datasheet for **TA384413**

### H3FT (HIST3H3) Mouse Monoclonal Antibody [Clone ID: 7G6-8F3-3H9]


#### Product data:

Product Type:	Primary Antibodies
Clone Name:	7G6-8F3-3H9
Applications:	WB
Recommended Dilution:	WB: 1/1000-3000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Synthetic Peptide of Histone H3 (Tri Methyl Lys36) (Methylated)
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.03% Proclin 300, pH 7.3.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Observed MW (kDa):15
Gene Name:	histone cluster 3, H3
Database Link:	<a href="#">Entrez Gene 8290 Human P68431</a>
Background:	Swiss-Prot Acc.P68431.Core component of nucleosome. Nucleosomes wrap and compact DNA into chromatin, limiting DNA accessibility to the cellular machineries which require DNA as a template. Histones thereby play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability. DNA accessibility is regulated via a complex set of post-translational modifications of histones, also called histone code, and nucleosome remodeling.MiscellaneousThis histone is only present in mammals and is enriched in acetylation of Lys-15 and dimethylation of Lys-10 (H3K9me2).
Synonyms:	H3.4; H3/g; H3/t; H3FT; H3t; MGC126886; MGC126888




[View online »](#)

Product images:



Western blot analysis of TriMethyl-Histone H3 (Lys36) (7G6) in THP-1, MCF-7, C6, HeLa lysates using TriMethyl-Histone H3 (Lys36) (7G6) antibody.



Western blot analysis of HeLa, diluted at 1) 1:1000 . 2) 1:5000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Invent biotech, MN, USA).