

## Product datasheet for **TA384420**

### H3FT (HIST3H3) Mouse Monoclonal Antibody [Clone ID: 5F6-6G6-6F9]

#### Product data:


Product Type:	Primary Antibodies
Clone Name:	5F6-6G6-6F9
Applications:	IF, IP, WB
Recommended Dilution:	WB: 1/1000-3000 ICC: 1/200 IP: 1/200
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Synthetic Peptide of Histone H3 (Di Methyl Lys27) (Methylated)
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.03% Proclin 300, pH 7.3.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Observed MW (kDa):15
Gene Name:	histone cluster 3, H3
Database Link:	<a href="#">Entrez Gene 8290 Human P68431</a>
Background:	Swiss-Prot Acc.P68431.Core component of nucleosome. Nucleosomes wrap and compact DNA into chromatin, limiting DNA accessibility to the cellular machineries which require DNA as a template. Histones thereby play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability. DNA accessibility is regulated via a complex set of post-translational modifications of histones, also called histone code, and nucleosome remodeling.MiscellaneousThis histone is only present in mammals and is enriched in acetylation of Lys-15 and dimethylation of Lys-10 (H3K9me2).




[View online »](#)

Synonyms: H3.4; H3/g; H3/t; H3FT; H3t; MGC126886; MGC126888


**Product images:**




Western blot detection of Di-Methyl-Histone H3 in C6,3T3,C2C12 cell lysates using Di-Methyl-Histone H3 (Lys27) (5F6) Mouse mAb(1:1000 diluted).Predicted band size:15kDa.Observed band size:15kDa.



Western blot detection of Di-Methyl-Histone H3 (Lys27) (5F6) in HeLa,Rat heart,Mouse heart,HUVEC,Jurkat cell lysates using Di-Methyl-Histone H3 (Lys27) (5F6) antibody. Observed band size:15kDa.



Immunoprecipitation analysis of DiMethyl-Histone H3 (Lys27) (5F6) in HeLa lysates using DiMethyl-Histone H3 (Lys27) (5F6) antibody.



Immunocytochemistry staining of HeLa cells using Di-Methyl-Histone H3 (Lys27) (5F6) Mouse mAb(dilution 1:100).