

## Product datasheet for TA384422S

### H3FT (HIST3H3) Mouse Monoclonal Antibody [Clone ID: 10F5-7A8-1B7]


#### Product data:

Product Type:	Primary Antibodies
Clone Name:	10F5-7A8-1B7
Applications:	WB
Recommended Dilution:	WB: 1/1000-3000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Synthetic Peptide of Histone H3 (Mono Methyl Lys9) (Methylated)
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.03% Proclin 300, pH 7.3.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Observed MW (kDa):15
Gene Name:	histone cluster 3, H3
Database Link:	<a href="#">Entrez Gene 8290 Human P68431</a>
Background:	Swiss-Prot Acc.P68431.Core component of nucleosome. Nucleosomes wrap and compact DNA into chromatin, limiting DNA accessibility to the cellular machineries which require DNA as a template. Histones thereby play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability. DNA accessibility is regulated via a complex set of post-translational modifications of histones, also called histone code, and nucleosome remodeling.MiscellaneousThis histone is only present in mammals and is enriched in acetylation of Lys-15 and dimethylation of Lys-10 (H3K9me2).
Synonyms:	H3.4; H3/g; H3/t; H3FT; H3t; MGC126886; MGC126888




[View online »](#)

Product images:



Western blot detection of Mono-Methyl-Histone H3 in MCF7,C2C12,C6,3T3 cell lysates using Mono-Methyl-Histone H3 (Lys9) (10F5) Mouse mAb(1:1000 diluted).Predicted band size:15KDa.Observed band size:15KDa.



Western blot detection of Mono-Methyl-Histone H3 (Lys9) in CHO-K1,Rat brain,C6,3T3,HeLa cell lysates using Mono-Methyl-Histone H3 (Lys9) (10F5) antibody. Observed band size:15kDa.