

## Product datasheet for **TA384445S**

### Heterogeneous Nuclear Ribonucleoprotein (A1 like) (HNRNPA1L2) Rabbit Monoclonal Antibody [Clone ID: R05-8G1]


#### Product data:

Product Type:	Primary Antibodies
Clone Name:	R05-8G1
Applications:	IHC, WB
Recommended Dilution:	WB: 1/1000-1/5000 IHC: 1/20-1/100
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A synthetic peptide of human Heterogeneous Nuclear Ribonucleoprotein (A1-like)
Formulation:	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Calculated MW: 34 kDa; Observed MW: 34 kDa
Gene Name:	heterogeneous nuclear ribonucleoprotein A1-like 2
Database Link:	<a href="#">Entrez Gene 144983 Human</a> <a href="#">Q5TBS2</a>
Background:	Swiss-Prot Acc.Q5TBS2.
Synonyms:	HNRNPA1L; LOC144983; MGC102957; OTTHUMP00000018460



[View online »](#)

## Product images:



Western blot analysis of Heterogeneous Nuclear Ribonucleoprotein (A1like) in K562, C6, 3T3, HeLa lysates using HNRNPA1L2 antibody.