


## Product datasheet for TA384450M

### hptII Mouse Monoclonal Antibody [Clone ID: 2D3-F3-D10]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	2D3-F3-D10
Applications:	WB
Recommended Dilution:	WB: 1/1000
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.03% Proclin 300, pH 7.3.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	37kDa

#### Product images:



Western blot analysis of HptII in Transgenic rice (leaf) and Nontransgenic rice (leaf) lysates using HptII antibody. Betaactin antibody (2000688F10) was used as a loading control.



[View online »](#)