

Product datasheet for TA384458

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

HSP70-1A (HSPA1A) Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: FC, IHC, IP, WB

Recommended Dilution: WB: 1/1000

IHC: 1/50 FC: 1/20 IP: 1/20

Reactivity: Human
Host: Rabbit
Isotype: IgG

Clonality: Polyclonal

Immunogen: Recombinant protein of human Hsp70

Formulation: 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA

Concentration: lot specific

Purification: Affinity Purified
Conjugation: Unconjugated

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Stability: 1 year

Predicted Protein Size: Calculated MW: 70 kDa; Observed MW: 70 kDa **Gene Name:** heat shock protein family A (Hsp70) member 1A

Database Link: Entrez Gene 3303 Human

P0DMV8

Background: Swiss-Prot Acc.P0DMV8.


Synonyms: FLJ54303; FLJ54370; FLJ54392; FLJ54408; FLJ75127; HSP70-1; HSP70-1/HSP70-2; HSP70-1A;

HSP70.1; HSP70I; HSP72; HSPA1; HSPA1B





Product images:



Immunohistochemistry of Hsp70 in paraffinembedded Human tonsil using Hsp70 Rabbit pAb at dilution 1/20