

Product datasheet for **TA384458S**


HSP70-1A (HSPA1A) Rabbit Polyclonal Antibody

Product data:

| | |
|--------------------------------|---|
| Product Type: | Primary Antibodies |
| Applications: | FC, IHC, IP, WB |
| Recommended Dilution: | WB: 1/1000 IHC: 1/50 FC: 1/20 IP: 1/20 |
| Reactivity: | Human |
| Host: | Rabbit |
| Isotype: | IgG |
| Clonality: | Polyclonal |
| Immunogen: | Recombinant protein of human Hsp70 |
| Formulation: | 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA |
| Concentration: | lot specific |
| Purification: | Affinity Purified |
| Conjugation: | Unconjugated |
| Storage: | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |
| Stability: | 1 year |
| Predicted Protein Size: | Calculated MW: 70 kDa; Observed MW: 70 kDa |
| Gene Name: | heat shock protein family A (Hsp70) member 1A |
| Database Link: | Entrez Gene 3303 Human P0DMV8 |
| Background: | Swiss-Prot Acc.P0DMV8. |
| Synonyms: | FLJ54303; FLJ54370; FLJ54392; FLJ54408; FLJ75127; HSP70-1; HSP70-1/HSP70-2; HSP70-1A; HSP70.1; HSP70I; HSP72; HSPA1; HSPA1B |



[View online »](#)

Product images:

Immunohistochemistry of Hsp70 in paraffin-embedded Human tonsil using Hsp70 Rabbit pAb at dilution 1/20