

Product datasheet for **TA384497**

IGHG1 Mouse Monoclonal Antibody [Clone ID: 2D7-A9-C8]


Product data:

Product Type:	Primary Antibodies
Clone Name:	2D7-A9-C8
Applications:	ELISA
Recommended Dilution:	ELISA: 1/10000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Purified recombinant human Human IgG protein fragments expressed in E.coli.
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.03% Proclin 300, pH 7.3.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year




[View online »](#)

Product images:



Recombinant Anthrax Protective Antigen used as Coating antigen, followed by AntiAnthrax Protective Antigen Human antibody and then incubated with AntiHuman IgG(201068) mouse mAb, HRP conjugated goat anti mouse IgG was used for coloring.



Recombinant Anthrax Protective Antigen used as Coating antigen, followed by different dilutions of supernatants which includes secreted AntiAnthrax Protective Antigen Human antibody, HRP conjugated AntiHuman IgG(201068) mouse mAb was used for coloring.