

Product datasheet for TA384520

IL1 beta (IL1B) Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: IF, IHC, WB

Recommended Dilution: WB: 1/500-1/2000

IHC: 1/50-1/200 IF: 1/50-1/200

Reactivity: Human, Mouse, Rat

Host: Rabbit Isotype: lgG

Clonality: Polyclonal

Immunogen: A synthetic peptide of human IL1B.Email for sequence

Formulation: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Concentration: lot specific

Purification: Affinity Purified Conjugation: Unconjugated

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. Storage:

Stability: 1 year

Predicted Protein Size: Calculated MW:30kDa;Observed MW:17kDa, 31kDa

Gene Name: interleukin 1 beta

Database Link: Entrez Gene 3553 Human

P01584

Background: Swiss-Prot Acc.P01584.The protein encoded by this gene is a member of the interleukin 1

> cytokine family. This cytokine is produced by activated macrophages as a proprotein, which is proteolytically processed to its active form by caspase 1 (CASP1/ICE). This cytokine is an important mediator of the inflammatory response, and is involved in a variety of cellular

> activities, including cell proliferation, differentiation, and apoptosis. The induction of cyclooxygenase-2 (PTGS2/COX2) by this cytokine in the central nervous system (CNS) is found to contribute to inflammatory pain hypersensitivity. This gene and eight other interleukin 1

family genes form a cytokine gene cluster on chromosome 2.



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn


Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com




Synonyms:

catabolin; IL-1; IL-1B; IL1-BETA; IL1F2; pro-interleukin-1-beta

Product images:





Immunofluorescence - IL1B Polyclonal Antibody. Immunofluorescence analysis of THP-1 cells using IL1B antibody at dilution of 1:100. Blue: DAPI for nuclear staining.




Immunofluorescence - IL1B Polyclonal Antibody. Immunofluorescence analysis of mouse lung using IL1B antibody at dilution of 1:100. Blue: DAPI for nuclear staining.









Western blot - IL1B Polyclonal Antibody. Western blot analysis of extracts of various cell lines, using IL1B antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Exposure time: 3s

Western blot - IL1B Polyclonal Antibody. Western blot analysis of extracts of THP-1 cells, using IL1B antibody at 1:1000 dilution. THP-1 cells were treated by PMA at 37°C for overnight. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25 up per lane. Blocking buffer: 3% nonfat dry milk in TBST. Exposure time: 3s.

Immunohistochemistry - IL1B Polyclonal Antibody. Immunohistochemistry of paraffinembedded rat lung using IL1B antibody at dilution of 1:1.