

## Product datasheet for **TA384529**

### IL7R alpha (IL7R) Rabbit Monoclonal Antibody [Clone ID: R08-5I2]


#### Product data:

Product Type:	Primary Antibodies
Clone Name:	R08-5I2
Applications:	WB
Recommended Dilution:	WB: 1/1000-1/5000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A synthetic peptide of human CD127
Formulation:	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Calculated MW: 52 kDa; Observed MW: 70 kDa
Gene Name:	interleukin 7 receptor
Database Link:	<a href="#">Entrez Gene 3575 Human P16871</a>
Background:	Swiss-Prot Acc.P16871.Receptor for interleukin-7. Also acts as a receptor for thymic stromal lymphopoietin (TSLP).
Synonyms:	CD127; CDW127; IL-7R-alpha; IL7RA; ILRA



[View online »](#)

## Product images:



Western blot analysis of CD127 in K562 lysates using CD127 antibody.