

## Product datasheet for **TA384586S**

### **JUNB Rabbit Monoclonal Antibody [Clone ID: R01-9B7]**


#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	R01-9B7
Applications:	IP, WB
Recommended Dilution:	WB: 1/1000 IP: 1/20
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A synthetic peptide of human JunB
Formulation:	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Calculated MW: 36 kDa; Observed MW: 42,43 kDa
Gene Name:	JunB proto-oncogene, AP-1 transcription factor subunit
Database Link:	<a href="#">Entrez Gene 3726 Human P17275</a>
Background:	Swiss-Prot Acc.P17275.Transcription factor involved in regulating gene activity following the primary growth factor response. Binds to the DNA sequence 5'-TGA[CG]TCA-3'.
Synonyms:	AP-1



[View online »](#)

## Product images:



Western blot analysis of JunB in K562, C6, 3T3, HeLa lysates using JunB antibody.