

Product datasheet for TA384693

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

MAGEA3 Rabbit Monoclonal Antibody [Clone ID: R07-8A2]

Product data:

Product Type: Primary Antibodies

Clone Name: R07-8A2

Applications: WB

Recommended Dilution: WB: 1/1000

Reactivity: Human
Host: Rabbit
Isotype: IgG

Clonality: Monoclonal

Immunogen: Recombinant protein of human MAGEA3

Formulation: 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA

Concentration: lot specific

Purification: Affinity Purified
Conjugation: Unconjugated

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Stability: 1 year

Predicted Protein Size: Calculated MW: 35 kDa; Observed MW: 45 kDa

Gene Name: MAGE family member A3

Database Link: Entrez Gene 4102 Human

P43357

Background: Swiss-Prot Acc.P43357.Proposed to enhance ubiquitin ligase activity of RING-type zinc finger-

containing E3 ubiquitin-protein ligases. May enhance ubiquitin ligase activity of TRIM28 and stimulate p53/TP53 ubiquitination by TRIM28. Proposed to act through recruitment and/or stabilization of the Ubl-conjugating enzyme (E2) at the E3:substrate complex. May play a role in embryonal development and tumor transformation or aspects of tumor progression. In vitro promotes cell viability in melanoma cell lines. Antigen recognized on a melanoma by


autologous cytolytic T-lymphocytes.

Synonyms: CT1.3; HIP8; HYPD; MAGE3; MAGEA6; MGC14613





Product images:



Western blot analysis of MAGEA3 in Hela lysates using MAGEA3 antibody.