

## Product datasheet for **TA384741M**


### **JNK1 (MAPK8) Rabbit Polyclonal Antibody**

#### **Product data:**

Product Type:	Primary Antibodies
Applications:	ELISA, IF, IHC, WB
Recommended Dilution:	IF: 1/50-200 WB: 1/500-2000
Reactivity:	Human, Mouse, Rat, Chicken
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized peptide derived from human JNK1/2/3 around the phosphorylation site of Thr183 and Tyr185. AA range:151-200 (Phosphorylated)
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Observed MW(kDa):46/54
Gene Name:	mitogen-activated protein kinase 8
Database Link:	<a href="#">Entrez Gene 5599 Human P45983</a>
Background:	Swiss-Prot Acc.P45983.
Synonyms:	JNK; JNK-46; JNK1; JNK1A2; JNK21B1/2; OTTHUMP00000019552; PRKM8; SAPK1



[View online »](#)

**Product images:**


Western blot analysis of Phospho-JNK (Thr183/Tyr185) in various lysates using Phospho-JNK (Thr183/Tyr185) antibody.