

## Product datasheet for **TA384776S**


### MGST3 Rabbit Monoclonal Antibody [Clone ID: R07-6A3]

#### Product data:


Product Type:	Primary Antibodies
Clone Name:	R07-6A3
Applications:	IHC, IP, WB
Recommended Dilution:	WB: 1/1000 IHC: 1/100 IP: 1/20
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A synthetic peptide of human MGST3
Formulation:	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Calculated MW: 17 kDa; Observed MW: 17 kDa
Gene Name:	microsomal glutathione S-transferase 3
Database Link:	<a href="#">Entrez Gene 4259 Human O14880</a>
Background:	Swiss-Prot Acc.O14880.Microsomal glutathione S-transferase 3
Synonyms:	GST-III; OTTHUMP00000032600



[View online »](#)

**Product images:**

Western blot analysis of MGST3 in K562, HeLa lysates using MGST3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using MGST3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.