

Product datasheet for TA384823

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

c-Myc (MYC) Rabbit Monoclonal Antibody [Clone ID: R05-5E2]

Product data:

Product Type: Primary Antibodies

Clone Name: R05-5E2

Applications: IF, IP, WB

Recommended Dilution: WB: 1/2000

ICC/IF: 1/20 IP: 1/20

Reactivity: Human, Mouse, Rat

Host: Rabbit Isotype: IgG

Clonality: Monoclonal

Immunogen: A synthetic peptide of human c-Myc

Formulation: 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA

Concentration: lot specific

Purification:Affinity PurifiedConjugation:Unconjugated

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Stability: 1 year

Predicted Protein Size: Calculated MW: 49 kDa; Observed MW: 57-65 kDa

Gene Name: v-myc avian myelocytomatosis viral oncogene homolog

Database Link: Entrez Gene 4609 Human


P01106

Background: Swiss-Prot Acc.P01106.
Synonyms: bHLHe39; c-Myc; MRTL






Product images:



Western blot analysis of c-Myc in K562, C6, 3T3, Hela lysates using c-Myc antibody.



Immunocytochemistry analysis of c-Myc (green) in Hela using c-Myc antibody,and DAPI(blue).