

## Product datasheet for **TA384834M**


### Visfatin (NAMPT) Rabbit Monoclonal Antibody [Clone ID: R05-4C4]

#### Product data:


Product Type:	Primary Antibodies
Clone Name:	R05-4C4
Applications:	IF, IHC, IP, WB
Recommended Dilution:	WB: 1/1000 IHC: 1/200 ICC/IF: 1/20 IP: 1/20
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A synthetic peptide of human Visfatin
Formulation:	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Calculated MW:56 kDa,Observed MW:56 kDa
Gene Name:	nicotinamide phosphoribosyltransferase
Database Link:	<a href="#">Entrez Gene 10135 Human P43490</a>
Synonyms:	1110035O14Rik; DKFZp666B131; MGC117256; NAMPTase; PBEF; PBEF1; VF; VISFATIN

[View online »](#)


## Product images:




Western blot detection of Visfatin in 3T3,Raw264.7 lysates using Visfatin antibody.Predicted band size:56kDa.Observed band size:56kDa.



Western blot detection of Visfatin in HeLa,A549,HL-60,U2OS,C6 cell lysates using Visfatin Rabbit mAb(1:1000 diluted).Predicted band size:56kDa.Observed band size:56kDa.



Immunohistochemistry of Visfatin in paraffin-embedded Human lung cancer tissue using Visfatin Rabbit mAb at dilution 1/50



Immunocytochemistry analysis of Nampt (green) in U87-MG using Nampt antibody, and DAPI (blue).