

## Product datasheet for **TA384888**


### **nptII Mouse Monoclonal Antibody [Clone ID: 4E9-C9-D11]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	4E9-C9-D11
Applications:	WB
Recommended Dilution:	WB: 1/1000
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.03% Proclin 300, pH 7.3.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	29kDa
Database Link:	<a href="#">P00552</a>



[View online »](#)

**Product images:**

Western blot analysis of nptII in nontransgenic Brassica napus L. and transgenic Brassica napus L. lysates using nptII antibody.