


Product datasheet for **TA384898**

NT5DC3 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1/500-2000
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthesized peptide derived from human NT5DC3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	60 kDa
Gene Name:	5'-nucleotidase domain containing 3
Database Link:	Entrez Gene 51559 Human Q86UY8
Background:	Swiss-Prot Acc.Q86UY8
Synonyms:	FLJ11266; GNN; TU12B1-TY; TU12B1TY

[View online »](#)

Product images:


Western blot analysis of NT5DC3 in A431 lysates using NT5DC3 antibody.