

## Product datasheet for **TA384905S**

### NUMB Rabbit Monoclonal Antibody [Clone ID: R09-2I5]

#### Product data:


Product Type:	Primary Antibodies
Clone Name:	R09-2I5
Applications:	IHC, IP, WB
Recommended Dilution:	WB: 1/1000 IHC: 1/200 IP: 1/20
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Recombinant protein of human NUMB
Formulation:	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Calculated MW: 71 kDa; Observed MW: 71 kDa
Gene Name:	NUMB, endocytic adaptor protein
Database Link:	<a href="#">Entrez Gene 8650 Human P49757</a>
Background:	Swiss-Prot Acc.P49757.Plays a role in the process of neurogenesis. Required throughout embryonic neurogenesis to maintain neural progenitor cells, also called radial glial cells (RGCs), by allowing their daughter cells to choose progenitor over neuronal cell fate. Not required for the proliferation of neural progenitor cells before the onset of neurogenesis. Also involved postnatally in the subventricular zone (SVZ) neurogenesis by regulating SVZ neuroblasts survival and ependymal wall integrity. May also mediate local repair of brain ventricular wall damage.




[View online »](#)

Synonyms: h-Numb; S171

### Product images:



Western blot analysis of NUMB in 3T3, C6, rat Brain, MCF-7 lysates using NUMB antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using NUMB antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.