

Product datasheet for **TA384910**

Omp Rabbit Monoclonal Antibody [Clone ID: R07-4I2]


Product data:

Product Type:	Primary Antibodies
Clone Name:	R07-4I2
Applications:	IP, WB
Recommended Dilution:	WB: 1/1000 IP: 1/20
Reactivity:	Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Recombinant protein of mouse Olfactory Marker Protein
Formulation:	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Calculated MW: 19 kDa; Observed MW: 19 kDa
Gene Name:	olfactory marker protein
Database Link:	Entrez Gene 18378 Mouse Q64288
Background:	Swiss-Prot Acc.Q64288.May act as a modulator of the olfactory signal-transduction cascade.
Synonyms:	OMP



[View online »](#)

Product images:



Western blot analysis of Olfactory Marker Protein in rat Brain, MCF-7 lysates using Olfactory Marker Protein antibody.