

Product datasheet for **TA384920S**


pagA Mouse Monoclonal Antibody [Clone ID: 5C1-C8-H8]

Product data:

Product Type:	Primary Antibodies
Clone Name:	5C1-C8-H8
Applications:	ELISA
Recommended Dilution:	ELISA: 1/10000
Reactivity:	Bacillus Anthracis
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.03% Proclin 300, pH 7.3.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	86kDa
Database Link:	P13423
Background:	Swiss-Prot Acc.P13423.



[View online »](#)

Product images:

Standard Curve for Anthrax Protective Antigen:
Capture Antibody mouse mAb 201240 to Anthrax
Protective Antigen at 4ug/ml and other antibody
for detecting.