

Product datasheet for **TA384944**


PDGF beta (PDGFB) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IF, IHC, WB
Recommended Dilution:	IF: 1/50-200 WB: 1/500-1/2000 IHC: 1/100-1/300 ELISA: 1/10000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized peptide derived from human PDGFB. AA range:16-65
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Observed MW (kDa):27
Gene Name:	platelet derived growth factor subunit B
Database Link:	Entrez Gene 5155 Human P01127
Background:	Swiss-Prot Acc.P01127.
Synonyms:	Becaplermin; c-sis; FLJ12858; PDGF-2; PDGF2; SIS; SSV



[View online »](#)

Product images:**Primary**
Antibody

— 抗

Primary Antibodies of Zen Bio were identified for Q&C tests by multiimmunology applications which were performed with different cell or tissue samples.