

## Product datasheet for **TA385017S**

### PPAR gamma (PPARG) Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IF, IHC, WB
Recommended Dilution:	IF: 1/50-200 WB: 1/500-1/2000 IHC: 1/100-1/300 ELISA: 1/10000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized peptide derived from human PPAR-gamma. AA range:78-127
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Observed MW (kDa):57
Gene Name:	peroxisome proliferator activated receptor gamma
Database Link:	<a href="#">Entrez Gene 5468 Human P37231</a>
Background:	Swiss-Prot Acc.P37231.
Synonyms:	CIMT1; GLM1; NR1C3; OTTHUMP00000185030; OTTHUMP00000185033; PPAR-gamma; PPARG1; PPARG2; PPARGgamma




[View online »](#)

## Product images:



Western blot analysis of PPAR gamma in HuvEc lysates using PPAR gamma antibody.



Immunohistochemical analysis of paraffin-embedded Rat-spleen tissue. PPAR- $\gamma$  Polyclonal Antibody was diluted at 1:200. Negative control was used by secondary antibody only.