

Product datasheet for TA385084

Ptgs2 Rabbit Monoclonal Antibody [Clone ID: R09-3C3]

Product data:

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn


Primary Antibodies
R09-3C3
IF, IP, WB
WB: 1/1000 ICC/IF: 1/20 IP: 1/20
Mouse
Rabbit
lgG
Monoclonal
Recombinant protein of mouse COX2
50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
lot specific
Affinity Purified
Unconjugated
Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
1 year
Calculated MW: 69 kDa; Observed MW: 74 kDa
prostaglandin-endoperoxide synthase 2
<u>Entrez Gene 19225 Mouse</u> <u>Q05769</u>
Swiss-Prot Acc.Q05769.Prostaglandin G/H synthase 2
COX-2; COX2; Cyclooxygenase-2; GRIPGHS; hCox-2; PGG/HS; PGHS-2; PHS-2



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US



Product images:



Western blot detection of COX2 in 3t3 cell lysates using COX2 Rabbit mAb(1:1000 diluted).Predicted band size:69KDa.Observed band size:74KDa.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US