

## Product datasheet for **TA385090M**


### FAK (PTK2) Rabbit Monoclonal Antibody [Clone ID: R05-1A1]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	R05-1A1
Applications:	WB
Recommended Dilution:	WB: 1/1000
Reactivity:	Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A synthetic phosphopeptide corresponding to residues surrounding Tyr397 of human FAK (Phosphorylated)
Formulation:	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Calculated MW: 119 kDa; Observed MW: 119 kDa
Gene Name:	protein tyrosine kinase 2
Database Link:	<a href="#">Q05397</a>
Background:	Swiss-Prot Acc.Q05397.
Synonyms:	FADK; FAK; FAK1; FRNK; pp125FAK

[View online »](#)

## Product images:



Western blot analysis of Phospho-FAK (Tyr397) in C6 lysates using Phospho-FAK (Tyr397) antibody.