

Product datasheet for TA385136

Rb (RB1) Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: ELISA, IF, IHC, WB

Recommended Dilution: WB: 1/500-1/2000

IHC: 1/100-1/300 IF: 1/200-1/1000 ELISA: 1/10000

Reactivity: Human, Mouse, Rat

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

Immunogen: The antiserum was produced against synthesized peptide derived from human

Retinoblastoma around the phosphorylation site of Ser780. AA range:751-800

(Phosphorylated)

Formulation: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Concentration: lot specific

Purification: Affinity Chromatography

Conjugation: Unconjugated

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Stability: 1 year

Predicted Protein Size: Observed MW (kDa):106

Gene Name: RB transcriptional corepressor 1

Database Link: Entrez Gene 5925 Human

P06400



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com




Background:

Swiss-Prot Acc.P06400.Key regulator of entry into cell division that acts as a tumor suppressor. Promotes G0-G1 transition when phosphorylated by CDK3/cyclin-C. Acts as a transcription repressor of E2F1 target genes. The underphosphorylated, active form of RB1 interacts with E2F1 and represses its transcription activity, leading to cell cycle arrest. Directly involved in heterochromatin formation by maintaining overall chromatin structure and, in particular, that of constitutive heterochromatin by stabilizing histone methylation. Recruits and targets histone methyltransferases SUV39H1, KMT5B and KMT5C, leading to epigenetic transcriptional repression. Controls histone H4 'Lys-20' trimethylation. Inhibits the intrinsic kinase activity of TAF1. Mediates transcriptional repression by SMARCA4/BRG1 by recruiting a histone deacetylase (HDAC) complex to the c-FOS promoter. In resting neurons, transcription of the c-FOS promoter is inhibited by BRG1-dependent recruitment of a phospho-RB1-HDAC1 repressor complex. Upon calcium influx, RB1 is dephosphorylated by calcineurin, which leads to release of the repressor complex .


Synonyms: OSRC; osteosarcoma; p105-Rb; pp110; pRb; RB

Product images:




Western blot analysis of Phospho-Rb (Ser780) in various lysates using Phospho-Rb (S780) antibody.






Western blot analysis of Phospho-Rb (Ser780) in KB lysates using Phospho-Rb (S780) antibody.



Western blot analysis of Phospho-Rb (Ser780) in K562 lysates treated with serum 10% using Phospho-Rb (Ser780) antibody. The lane on the right is blocked with the Phospho-peptide.



Immunohistochemistry analysis of paraffinembedded Human breast carcinoma using Retinoblastoma (Phospho-Ser78) antibody. Highpressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Sample with blocking peptide on the right.