

Product datasheet for **TA385175S**

RPA70 (RPA1) Mouse Monoclonal Antibody [Clone ID: 8C3-D12-H10]


Product data:

Product Type:	Primary Antibodies
Clone Name:	8C3-D12-H10
Applications:	IF, IP, WB
Recommended Dilution:	WB: 1/1000 ICC: 1/100
Reactivity:	Human, Monkey, Mouse, Rat
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Purified recombinant human RPA70 protein fragments expressed in E.coli.
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.03% Proclin 300, pH 7.3.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	70kDa
Gene Name:	replication protein A1
Database Link:	Entrez Gene 6117 Human P27694
Background:	Swiss-Prot Acc.P27694.
Synonyms:	HSSB; MST075; MSTP075; REPA1; RF-A; RP-A; RPA70



[View online »](#)

Product images:



Western blot detection of RPA70 in C6, COS7, K562 and HeLa cell lysates using RPA70 mouse mAb (1:1000 diluted). Predicted band size: 70 kDa. Observed band size: 70 kDa.