

Product datasheet for **TA385183**


RSK1 p90 (RPS6KA1) Mouse Monoclonal Antibody [Clone ID: 6B9-D6-F8]

Product data:


Product Type:	Primary Antibodies
Clone Name:	6B9-D6-F8
Applications:	IP, WB
Recommended Dilution:	WB: 1/1000 IP: 1/20
Reactivity:	Human, Mouse, Monkey
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Purified recombinant human p90RSK protein fragments expressed in E.coli.
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.03% Proclin 300, pH 7.3.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	90kDa
Gene Name:	ribosomal protein S6 kinase A1
Database Link:	Entrez Gene 6195 Human Q15418
Background:	Swiss-Prot Acc.Q15418.Serine/threonine kinase that may play a role in mediating the growth-factor and stress induced activation of the transcription factor CREB.
Synonyms:	HU-1; MAPKAPK1A; p90S6K; pp90RSK1; RSK; RSK-1; RSK1



[View online »](#)

Product images:


Immunoprecipitation analysis of HeLa cell lysates using p90RSK mouse mAb.



Western blot analysis of p90RSK in COS7, 3T3, Jurkat, HeLa and A549 lysates using p90RSK antibody.