

## Product datasheet for **TA385297S**

### SH3PX1 (SNX9) Rabbit Monoclonal Antibody [Clone ID: R03-5I5]


#### Product data:

Product Type:	Primary Antibodies
Clone Name:	R03-5I5
Applications:	IF, IHC, IP, WB
Recommended Dilution:	WB: 1/1000 IHC: 1/50 ICC/IF: 1/20 IP: 1/20
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Recombinant protein of human SH3PX1
Formulation:	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Calculated MW: 67 kDa; Observed MW: 67 kDa
Gene Name:	sorting nexin 9
Database Link:	<a href="#">Entrez Gene 51429 Human Q9Y5X1</a>
Background:	Swiss-Prot Acc.Q9Y5X1.Sorting nexin-9
Synonyms:	MST155; MSTP155; SDP1; SH3PX1; SH3PXD3A; WISP




[View online »](#)


Product images:




Western blot analysis of SNX9 in mouse intestines lysates using SNX9 antibody.



Western blot analysis of SH3PX1 in K562, HeLa lysates using SH3PX1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using SNX9 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of SNX9 (green) in HeLa using SNX9 antibody, and DAPI (blue).