

## **Product datasheet for TA385301**

### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## Superoxide Dismutase 1 (SOD1) Rabbit Monoclonal Antibody [Clone ID: R07-1G6]

#### **Product data:**

**Product Type:** Primary Antibodies

**Clone Name:** R07-1G6 **Applications:** IHC, WB

Recommended Dilution:

WB: 1/1000 IHC: 1/20

Reactivity: Human, Mouse, Rat

Host: Rabbit Isotype: IgG

Clonality: Monoclonal

Immunogen: A synthetic peptide of human Superoxide Dismutase 1

Formulation: 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA

**Concentration:** lot specific

Purification: Affinity Purified Conjugation: Unconjugated

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Stability: 1 year

Predicted Protein Size: Calculated MW: 16 kDa; Observed MW: 18 kDa

**Gene Name:** superoxide dismutase 1, soluble

Database Link: Entrez Gene 6647 Human

P00441

Background: Swiss-Prot Acc.P00441.Destroys radicals which are normally produced within the cells and

which are toxic to biological systems. Miscellaneous The protein (both wild-type and ALS1 variants) has a tendency to form fibrillar aggregates in the absence of the intramolecular disulfide bond or of bound zinc ions. These aggregates may have cytotoxic effects. Zinc


binding promotes dimerization and stabilizes the native form.

Synonyms: ALS1; homodimer; IPOA; SOD






# **Product images:**



Western blot analysis of Superoxide Dismutase 1 in K562, rat Brain, C6, 3T3, Hela lysates using Superoxide Dismutase 1 antibody.



Immunohistochemistry analysis of paraffinembedded Human colon cancer using Superoxide Dismutase 1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.