

## **Product datasheet for TA385306**

## 1 Todact datasticct for TASOSSO

## Product data:

**Product Type:** Primary Antibodies

**Applications:** IF, IHC, WB

**SP1 Rabbit Polyclonal Antibody** 

Recommended Dilution: WB: 1/500-1/2000

IHC: 1/50-1/200 ICC/IF: 1/50-1/200

Reactivity: Human
Host: Rabbit
Isotype: IgG

Clonality: Polyclonal

**Immunogen:** A synthesized peptide derived from human SP1

Formulation: Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50%

glycerol.

**Concentration:** lot specific

**Purification:** Affinity Chromatography

**Conjugation:** Unconjugated

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Stability: 1 year Predicted Protein Size: 90kDa

Gene Name: Sp1 transcription factor

Database Link: Entrez Gene 6667 Human

P08047

**Background:** Swiss-Prot Acc.P08047.SP1 is a transcription factor of the Sp1 C2H2-type zinc-finger protein

family. Phosphorylated and activated by MAPK. Dephosphorylation by PTEN inhibits DNA

binding. Binds to p38 in the nucleus.

Synonyms: TSFP1




**OriGene Technologies, Inc.** 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



## **Product images:**



Western blot analysis of SP1 in HeLa lysates using SP1 antibody.