

Product datasheet for **TA385315M**

SQSTM1 Rabbit Monoclonal Antibody [Clone ID: R08-4H8]


Product data:

Product Type:	Primary Antibodies
Clone Name:	R08-4H8
Applications:	IF, IP, WB
Recommended Dilution:	WB: 1/1000 ICC/IF: 1/20 IP: 1/20
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Recombinant protein of human SQSTM1/p62
Formulation:	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Calculated MW:48 kDa,Observed MW:62 kDa
Gene Name:	sequestosome 1
Database Link:	Entrez Gene 8878 Human Q13501
Synonyms:	A170; EBIAP; ORCA; OSIL; p60; p62; p62B; PDB3; ZIP3




[View online »](#)

Product images:



Western blot detection of SQSTM1/p62 in A549,HL-60,U2OS cell lysates using SQSTM1/p62 Rabbit mAb(1:1000 diluted).Predicted band size:48kDa.Observed band size:62kDa.



Immunocytochemistry analysis of SQSTM1/p62 (green) in U87-MG using SQSTM1/p62 antibody,and DAPI(blue).