

Product datasheet for **TA385332**

STAT2 Rabbit Polyclonal Antibody


Product data:

| | |
|-------------------------|--|
| Product Type: | Primary Antibodies |
| Applications: | IF, IHC, WB |
| Recommended Dilution: | WB: 1/2000 IHC: 1/50 ICC/IF: 1/20 |
| Reactivity: | Human |
| Host: | Rabbit |
| Isotype: | IgG |
| Clonality: | Polyclonal |
| Immunogen: | A synthetic peptide of human STAT2 |
| Formulation: | 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA |
| Concentration: | lot specific |
| Purification: | Affinity Purified |
| Conjugation: | Unconjugated |
| Storage: | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |
| Stability: | 1 year |
| Predicted Protein Size: | Calculated MW: 98 kDa; Observed MW: 98 kDa |
| Gene Name: | signal transducer and activator of transcription 2 |
| Database Link: | Entrez Gene 6773 Human P52630 |
| Background: | Swiss-Prot Acc.P52630. |
| Synonyms: | ISGF-3; MGC59816; p113; STAT113 |



[View online »](#)

Product images:



Western blot analysis of STAT2 in K562 lysates using STAT2 antibody.