

Product datasheet for **TA385351**

Stathmin 1 (STMN1) Mouse Monoclonal Antibody [Clone ID: 10D6-2D10-2F9]

Product data:

Product Type:	Primary Antibodies
Clone Name:	10D6-2D10-2F9
Applications:	IHC
Recommended Dilution:	IHC-p: 1/100-1/500 WB: 1/500-1/2000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2b, kappa
Clonality:	Monoclonal
Immunogen:	Synthesized peptide derived from human Stathmin
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Gene Name:	stathmin 1
Database Link:	Entrez Gene 3925 Human P16949
Background:	Swiss-Prot Acc.P16949.
Synonyms:	C1orf215; FLJ32206; Lag; LAP18; metablastin; MGC138869; MGC138870; OP18; OTTHUMP00000008530; PP17; PP19; PR22; prosolin; SMN



[View online »](#)

Product images:

