

## Product datasheet for **TA385356S**

### Striatin 4 (STRN4) Rabbit Monoclonal Antibody [Clone ID: R04-4J1]


#### Product data:

Product Type:	Primary Antibodies
Clone Name:	R04-4J1
Applications:	IF, IP, WB
Recommended Dilution:	WB: 1/1000-1/5000 ICC/IF: 1/20 IP: 1/20
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Recombinant protein of human Striatin 4
Formulation:	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Calculated MW: 81 kDa; Observed MW: 81 kDa
Gene Name:	striatin 4
Database Link:	<a href="#">Entrez Gene 29888 Human Q9NRL3</a>
Background:	Swiss-Prot Acc.Q9NRL3.Binds calmodulin in a calcium dependent manner. May function as scaffolding or signaling protein.MiscellaneousThe name 'Zinedin' probably originates from the name of the famous soccer player from Marseille (Zinedine Zidane).
Synonyms:	FLJ35594; ZIN; zinedin



[View online »](#)

## Product images:



Western blot analysis of Striatin 4 in HeLa lysates using Striatin 4 antibody.