

### **Product datasheet for TA385401**

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## TGF beta Receptor II (TGFBR2) Rabbit Polyclonal Antibody

#### **Product data:**

Product Type: Primary Antibodies

Applications: ELISA, IF, IHC, WB

Recommended Dilution: IF: 1/50-200

WB: 1/500-1/2000 IHC: 1/100-1/300 ELISA: 1/20000

Reactivity: Human, Mouse, Rat

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

Immunogen: The antiserum was produced against synthesized peptide derived from human TGF beta

Receptor II. AA range:91-140

**Formulation:** Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

**Concentration:** lot specific

Purification: Affinity Purified
Conjugation: Unconjugated

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Stability: 1 year

Predicted Protein Size: Observed MW (kDa):65

**Gene Name:** transforming growth factor beta receptor 2

Database Link: Entrez Gene 7048 Human


P37173

**Background:** Swiss-Prot Acc.P37173.

Synonyms: AAT3; FAA3; HNPCC6; LDS1B; LDS2B; MFS2; RIIC; TAAD2; TbetaR-II; TGFbeta-RII; TGFR-2



# **Product images:**



Western blot analysis of TGF beta Receptor II in Jurkat lysates using TGF beta Receptor II antibody.