

Product datasheet for **TA385413**


Tlr13 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 1/200-500
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Protein of TLR13
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Gene Name:	toll-like receptor 13
Database Link:	Entrez Gene 279572 Mouse Q6R5N8
Background:	Swiss-Prot Acc.Q6R5N8.
Synonyms:	AI666735; Gm713



[View online »](#)

Product images:

Immunohistochemistry analysis of paraffin-embedded mouse Spleen Tissue using Toll-Like Receptor 13 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.