

Product datasheet for **TA385415M**

TLR3 Mouse Monoclonal Antibody [Clone ID: 2F10-B8-G9]


Product data:

Product Type:	Primary Antibodies
Clone Name:	2F10-B8-G9
Applications:	WB
Recommended Dilution:	WB: 1/1000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	Purified recombinant human Toll-Like Receptor 3 protein fragments expressed in E.coli.
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.03% Proclin 300, pH 7.3.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Refer to figures
Gene Name:	toll like receptor 3
Database Link:	Entrez Gene 7098 Human O15455
Background:	Swiss-Prot Acc.O15455.
Synonyms:	CD283



[View online »](#)

Product images:



Western blot analysis of Toll-Like Receptor 3 in CHO-K1 lysates and CHO-K1 transfected by Toll-Like Receptor 3 lysates using Toll-Like Receptor 3 antibody.