

Product datasheet for **TA385423M**


TMEM9B Rabbit Monoclonal Antibody [Clone ID: R02-9C6]

Product data:


Product Type:	Primary Antibodies
Clone Name:	R02-9C6
Applications:	WB
Recommended Dilution:	WB: 1/1000-1/2000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Recombinant protein of human TMEM9B
Formulation:	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Calculated MW: 23 kDa; Observed MW: 23 kDa
Gene Name:	TMEM9 domain family member B
Database Link:	Entrez Gene 56674 Human Q9NQ34
Background:	Swiss-Prot Acc.Q9NQ34.Enhances production of proinflammatory cytokines induced by TNF, IL1B, and TLR ligands. Has a role in TNF activation of both the NF-kappaB and MAPK pathways.
Synonyms:	C11orf15

[View online »](#)

Product images:



Western blot analysis of TMEM9B in Hacat, PC-12 lysates using Transmembrane Protein 9B antibody.



Western blot analysis of TMEM9B in 3T3 lysates using TMEM9B antibody.