

Product datasheet for **TA385426S**

TNF alpha (TNF) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IF, IHC, WB
Recommended Dilution:	WB: 1/500-1/2000 IHC: 1/100-1/300 IF: 1/200-1/1000 ELISA: 1/10000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized peptide derived from human TNFA. AA range:141-190
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Observed MW (kDa):15 - 30
Gene Name:	tumor necrosis factor
Database Link:	Entrez Gene 7124 Human P01375
Background:	Swiss-Prot Acc.P01375.
Synonyms:	Cachectin; DIF; TNF-a; TNF-alpha; TNFA; TNFSF2



[View online »](#)

Product images:



Western Blot analysis of Jurkat cells using TNF-u03b1 Polyclonal Antibody