

Product datasheet for **TA385432S**

TNFRSF1B Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	FC, IF, IHC, IP, WB
Recommended Dilution:	WB: 1/2000 IHC: 1/50 ICC/IF: 1/50 FC: 1/20 IP: 1/20
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	A synthetic peptide of human TNF Receptor II
Formulation:	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Calculated MW: 48 kDa; Observed MW: 65 kDa
Gene Name:	tumor necrosis factor receptor superfamily member 1B
Database Link:	Entrez Gene 7133 Human P20333
Background:	Swiss-Prot Acc.P20333.
Synonyms:	CD120b; Etanercept; p75; p75TNFR; TBPII; TNF-R-II; TNF-R2; TNF-R75; TNF-RII; TNFBR; TNFR-II; TNFR1B; TNFR2; TNFR80




[View online »](#)


Product images:




Western blot analysis of TNF Receptor II in 3T3, C6 lysates using TNF Receptor 2 antibody.



Western blot analysis of TNF Receptor II in Jurkat, HeLa lysates using TNF Receptor II antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using TNF Receptor II antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of TNF Receptor II (green) in HeLa using TNF Receptor II antibody, and DAPI (blue)