

Product datasheet for **TA385443M**

p53 (TP53) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	WB: 1/500-1/2000 IHC-p: 1/100-300 ELISA: 1/20000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthesized acetyl-peptide derived from the human p53 around the acetylation site of K382. (Acetylated)
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Observed MW (kDa):53
Gene Name:	tumor protein p53
Database Link:	Entrez Gene 7157 Human P04637
Background:	Swiss-Prot Acc.P04637.
Synonyms:	FLJ92943; LFS1; p53; TRP53



[View online »](#)

Product images:



Western blot analysis of Acetyl-p53 (Lys382) in HepG2 lysates using Acetyl-p53 (Lys382) antibody.