

Product datasheet for **TA385454**

TCTP (TPT1) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	FC, IF, IHC, IP, WB
Recommended Dilution:	WB: 1/1000-1/5000 IHC: 1/50-1/200 ICC/IF: 1/50-1/200 IP: 1/50 FC: 1/100
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	A synthesized peptide derived from human TCTP
Formulation:	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	20kDa
Gene Name:	tumor protein, translationally-controlled 1
Database Link:	Entrez Gene 7178 Human P13693
Background:	Swiss-Prot Acc.P13693.Involved in calcium binding and microtubule stabilization.
Synonyms:	FLJ27337; Fortilin; HRF; p02; p23; TCTP



[View online »](#)

Product images:



Western blot analysis of TCTP Receptor in HeLa lysates using TCTP antibody.