

## Product datasheet for **TA385509**

### Ubiquitin (UBB) Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	WB: 1/500-1/2000 IHC: 1/100-1/300 ELISA: 1/20000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized peptide derived from human Ubiquitin. AA range:40-89
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Refer to figures
Gene Name:	ubiquitin B
Database Link:	<a href="#">Entrez Gene 7314 Human P0CG47</a>
Background:	Swiss-Prot Acc.P0CG47/P0CG48/P62979/P62987.
Synonyms:	FLJ25987; MGC8385; OTTHUMP00000064960; OTTHUMP00000064961; UBA80; UBCEP1; UBCEP2



[View online »](#)

## Product images:



Western blot analysis of Ubiquitin in various lysates using Ubiquitin antibody.