

Product datasheet for **TA385521**

UFC1 Rabbit Monoclonal Antibody [Clone ID: R08-5E1]


Product data:

Product Type:	Primary Antibodies
Clone Name:	R08-5E1
Applications:	IHC, WB
Recommended Dilution:	WB: 1/2000-1/10000 IHC: 1/20
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Recombinant protein of human UFC1
Formulation:	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Calculated MW: 20 kDa; Observed MW: 20 kDa
Gene Name:	ubiquitin-fold modifier conjugating enzyme 1
Database Link:	Entrez Gene 51506 Human Q9Y3C8
Background:	Swiss-Prot Acc.Q9Y3C8.E2-like enzyme which forms an intermediate with UFM1 via a thioester linkage.
Synonyms:	HSPC155




[View online »](#)

Product images:



Western blot analysis of UFC1 in K562, rat Brain, C6, 3T3 lysates using UFC1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using UFC1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.