

## Product datasheet for **TA385531**

### USP11 Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	FC, IF, IP, WB
Recommended Dilution:	WB: 1/500-1/2000 ICC/IF: 1/50-1/200 IP: 1/100 FC: 1/200
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	A synthesized peptide derived from human USP11
Formulation:	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	110kDa
Gene Name:	ubiquitin specific peptidase 11
Database Link:	<a href="#">Entrez Gene 8237 Human P51784</a>
Background:	Swiss-Prot Acc.P51784. Protease that can remove conjugated ubiquitin from target proteins and polyubiquitin chains. Inhibits the degradation of target proteins by the proteasome. Plays a role in the regulation of pathways leading to NF-kappa-B activation. Plays a role in the regulation of DNA repair after double-stranded DNA breaks.
Synonyms:	UHX1



[View online »](#)

## Product images:



Western blot analysis of USP11 in (1) HEK293 lysates; (2) mouse testis lysates using USP11 antibody.