

## Product datasheet for **TA385595**

### Insulin Glargine Mouse Monoclonal Antibody [Clone ID: 7A6-6G6-9G5]

#### Product data:

|                       |  |
|-----------------------|--|
| Product Type:         | Primary Antibodies   |
| Clone Name:           | 7A6-6G6-9G5  |
| Applications:         | ELISA  |
| Recommended Dilution: | ELISA: 1/2000-10000  |
| Host:                 | Mouse  |
| Isotype:              | IgG1   |
| Clonality:            | Monoclonal   |
| Immunogen:            | Insulin Glargine   |
| Formulation:          | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.03% Proclin 300, pH 7.3.           |
| Concentration:        | lot specific   |
| Purification:         | Affinity Purified  |
| Conjugation:          | Unconjugated   |
| Storage:              | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |
| Stability:            | 1 year   |



[View online »](#)