

## Product datasheet for **TA385597**

### Rubisco Mouse Monoclonal Antibody [Clone ID: 1E3-6F10-5H9]


#### Product data:

Product Type:	Primary Antibodies
Clone Name:	1E3-6F10-5H9
Applications:	WB
Recommended Dilution:	WB: 1/2000-10000
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Purified recombinant protein expressed in E.coli.
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.03% Proclin 300, pH 7.3.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	53kDa



[View online »](#)

## Product images:



Western blot analysis of Rubisco (1E3) in Arabidopsis lysates using Rubisco (1E3) antibody.